P0.48

**Heat Pump,
Compressor Principle
04370.88**

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
Fig. 1

1. PURPOSE AND CHARACTERISTIC FEATURES

The heat pump is the demonstration model of a compressor heat pump. Pressure and temperature of the refrigerant can be measured on it and its state observed in sight glasses.

2. PRINCIPLES

2.1 The concept of a heat pump

The purpose of a heat pump is to extract the quantity of heat Q_1 from an environment at a low temperature T_1 and to deliver it to another environment having a higher temperature T_2 . Work W is necessary for this. This work is carried out by a compressor. A compressor heat pump works on the principle of a cyclic process, which should approximate as far as possible to an ideal Carnot cycle. See in this respect the TS diagram in Fig. 2; the absolute temperature T is recorded over the entropy $S = Q/T$.

In detail, the cyclic process operates as follows:

Starting at point 1, the vapourised working medium at pressure p_1 and temperature T_1 is compressed at point 2 to pressure p_2 , where-

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
Fig. 3

line 2 and temperature measuring point 3 into the condenser 4. There, it releases the surplus heat to the environment, e.g. water in a vessel (not shown). The refrigerant thereby condenses and passes through the temperature measuring point 5 into the red line 6. The aggregate state of the refrigerant prevailing in this line can be observed in the sight glass 7. The refrigerant pressure is read off on the pressure gauge 8. Through line 9 the refrigerant reaches the throttle valve 13, is pressure-relieved there and thereby cools. It flows through the blue line 14 and temperature measuring point 15 to the evaporator 16. There, it takes up heat from the environment, becomes gaseous and flows through temperature measuring point 17 into the dotted blue line 18. The refrigerant now has a heat content Q_k . The gaseous state of the refrigerant can be seen in the sight glass 11. The pressure gauge 12 indicates the pressure of the gaseous refrigerant. This now returns via line 10 to the compressor 1, and the process starts anew.

Some components, which become visible when the rear wall (fastened with 4 knurled screws) is removed, are of importance as regards satisfactory operation of the heat pump but not for explaining its method of operation:

- A receiver with a store of refrigerant to compensate for losses by possible leaks
- A drier which removes any water residues from the refrigerant
- A high-low pressure switch to protect the heat pump against overheating or undercooling; the pressure switch temporarily switches off the compressor outside certain pressure limits set in the factory.

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Fig. 2: Diagramme T, S du cycle Carnot

Fig. 2: Diagrama T-S del ciclo de Carnot.

by the temperature of the medium is increased from T_1 to T_2 . This involves no exchange of heat energy with the environment, so that $dQ=0$ and $S=\text{const.}$

At point 2, the outlet valve of the compressor opens and the working medium flows into the condenser. There, it becomes liquid and gives off, at temperature T_2 , the quantity of heat Q_2 as useful heat. The temperature T_2 remains unchanged until condensation is complete. Because $dS = dQ/T$ and $T = \text{const.}$, the entropy S steadily decreases until it has fallen from S_b to S_a (at point 3 in the diagram).

From point 3 to point 4 the liquid medium flows through the throttling device and its pressure is thereby reduced to p_1 . The temperature also drops from T_2 to T_1 , the entropy remaining unchanged.

The working medium flows into the evaporator without absorbing or releasing heat. It there absorbs the quantity of heat Q_1 from the heat source at temperature T_1 and vapourises. The vapour is drawn off by the compressor. The energy stored in the working medium increases until the medium is fully vapourised. As the temperature does not change during this process, the entropy increases from S_a to S_b . The process has now reached the starting point at 1 and the procedure recommences.

The area 1-2-3-4 in the diagram represents the work W performed by the compressor, and the area 1-4- S_a - S_b represents the heat extracted from the heat source Q_1 . In the ideal case portrayed above, the area S_b -2-3- S_a is then the useful heat Q_2 obtained at temperature T_2 , as $Q_2 = Q_1 + W$.

3. DESCRIPTION AND MODE OF OPERATION

The heat pump is designed as a demonstration unit in which various component parts are integrated into a schematic diagram in order to give a clear idea of the cyclic process taking place (see Fig. 3). The closed system contains refrigerant R 12 as the working medium. The compressor 1, the pipelines for which run hidden behind the supports of the stand, compresses the gaseous refrigerant (coming from the evaporator 16) and increases the available quantity of heat Q_k by Q_e to Q_w . The superheated, gaseous refrigerant is forced through the dotted red

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4. OPERATING INSTRUCTIONS

Connect the connecting lead of the heat pump to the AC mains; the mains switch is on the side of the casing.

Warning! On delivery of the unit, to protect the compressor wait about an hour before first starting up, as the distribution of refrigerant in the system can be disturbed during transport by shaking or by the unit not being in a vertical position.

Immerse both heat exchangers in water containers in order to carry out the tests. To do this, the unit should either be raised (avoiding undue tipping) or moved so that the heat exchangers protrude one after the other over the edge of the table, and the containers (which are best filled only then) can be pushed underneath.

The temperature probes intended (see Equipment List) should be set into the temperature measuring points, using a little heat-conductive paste; additional fastening is not necessary.

The heat exchangers should be freed of lime deposits from time to time; diluted acetic acid (about 5%) is suitable for this.

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6. TECHNICAL DATA

Standard rating of compressor 150 W
Efficiency $\eta = 81\%$
Coefficient of performance $e = 2.6$
(maximum)
Operating voltage 220 V; 50 Hz

7. EQUIPMENT LIST

04370.88 Heat pump, compressor principle, including 2 insulating containers with cock; capacity 5 l

For temperature measurement optionally (4 measuring probes required in each case)

38005.00 Thermometer $-10^\circ \dots +100^\circ\text{C}$

or for demonstration measurement

11100.00 Moving coil instrument

11106.02 Range multiplier $+120^\circ\text{C}$

07922.01 Round cell, 1.5 V, R 14 (2 required)

04187.00 Temperature probe

or

11759.93 Digital apparatus for temperature measurement 4-4

11759.01 Temperature probe, immersion type, Pt 100

or

11762.93 Digital apparatus for temperature measurement 4-1

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PHYWE AG