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1. Description

The double wire loop oscillograph is used to examine low-frequency electric oscillations, in particular the alternating current of 50 hertz. Since the apparatus has two loops which are independent of each other, it is possible to trace two different curves at the same time and to examine them for their phase difference.

The natural frequency of the loops is approximately 70 hertz; the sensitivity for readings with spot of light is 3—5 mA/mm/m. The loops can be loaded up to 300 mA although in many cases much smaller currents suffice.

Fig. 1 gives a schematic illustration of the apparatus. The loops are stretched between pole shoes (P). They can be connected up independent of each other, each by means of a pair of terminals (K); each loop has a concave mirror (H). By means of the adjusting device (V), it is pos-

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
Notes

1. The five-figure numbers quoted in brackets refer to the catalogue numbers of the respective apparatus.
2. The data DK . . . refer to the experiment descriptions published as a collection of leaflets in the "LEYBOLD PHYSIKALISCHE HANDBLÄTTER".
3. The data and illustrations are not binding in every detail for the design of the apparatuses. It is our sole aim always to adapt our manufacturing programme to the most recent knowledge gained in all scientific and technical fields.

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sible to change the position of the loops, the mirror and the images (see below).

To carry out the experiment, the apparatus should be mounted on a U-core (562 11) which produces the necessary magnetic field.

2. Assembly and illumination

Place two coils of 250 turns (562 13) on a U-core (562 11). Remove the clamping device, and in its place mount the double wire loop oscillograph by screwing it on tightly so that the pole shoes lie on the surfaces of the U-core arms.

The assembled apparatus should be set up with the loops horizontal, and illuminated so that the light ray emitted by the lamp is reflected once from the concave mirror and once from a rotating mirror cube (415 21) or from the rotating mirror (504 96) respectively, falling onto a screen or a smooth wall. The mirror (504 96) should be mounted on the large electric stop-clock, and the mirror cube (415 21) on stand equipment. The light is best produced by means of the lamp housing (450 60 a) with lamp (450 51), and the single-lens condenser with diaphragm (460 17). The lamp coil should be set in the horizontal, and projected onto the concave mirrors by the condenser. Its distance from the double wire loop oscillograph should be selected so that the concave mirrors produce well-defined images of the aperture in the diaphragm. These images should then be illuminated evenly by turning the adjusting screws on the lamp housing. The two images can easily be brought to coincidence by manipulating the adjusting device (V). The level should be corrected by turning, the lateral deviations by sliding the device.

The complete illuminating apparatus is supplied under Cat. No. 516 02. It consists of the following parts: 1 stand base, 20 cm (300 02), 1 stand column, 50 cm long (300 42), 1 LEYBOLD's universal clamp (301 01), lamp housing (450 60 a) with lamp (450 51) and single-lens condenser (460 17).

3. Electric connection

The magnetic field should be excited with D. C. taken from a rectifier (e. g. 523 64), the switchboard or an accumulator. Approximately 5 A are required when using coils with 250 turns, and 2.5 A for coils with 500 turns. To pole the coils correctly, connect up the sockets marked (E), and

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4. Versuche

in DK-Abt. 537 und 538.

connect the sockets (A) to the source of current.

The wiring illustrated in Fig. 2 is suited for demonstrating phase difference. The voltage of approximately 50 volts should be taken from an experiment transformer (562 11/12) with the coils with 1,000 and 250 turns. By using the sliding rheostat, the size of the two currents can be varied, so that one can differentiate between the curves which appear when the mirror is rotated.

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