
**P67.16**

**Faraday-Becher**  
**Faraday's Ice Pail**  
**Cylindre de Faraday**

1:4

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**1. Description**

The Faraday's ice pail is used to examine the distribution of charge on a tubular conductor. Moreover it shows that the inside of a charged conductor always remains uncharged and that by feeding many small quantities of charge the pail is charged to high voltages.

The pail is made of aluminium. It is about 13 cm long with a diameter of approximately 7.5 cm. A plug affixed to the base enables it to be mounted either on the simple electroscop (540 09) or on the Wulf's electroscop (546 01 a).

**2. Usage**

The Faraday's ice pail can be employed in several ways:

- 1) Place the pail on an insulated base e. g. paraffin or Holtz's terminal insulator (540 50), and charge it by means of the Van de Graaf generator (541 69) or the Wimshurst machine (541 63).

With the aid of a "charge carrier" (metal disk 542 52 or test condenser 546 23) and an electroscop it can be demonstrated that no charge can be taken away from

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**Notes**

1. The five-figure numbers quoted in brackets refer to the catalogue numbers of the respective apparatus.
2. The data DK . . . refer to the experiment descriptions published as a collection of leaflets in the "LEYBOLD PHYSIKALISCHE HANDBLÄTTER".
3. The data and illustrations are not binding in every detail for the design of the apparatuses. It is our sole aim always to adapt our manufacturing programme to the most recent knowledge gained in all scientific and technical fields.

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the inside of the charged pail. In carrying out this demonstration touch the base only and not the sides of the pail because inside its upper part are still several lines of force getting in from outside. It is particularly impressive to mount the pail direct on the Van de Graaf generator.

- 2) Place the pail on top of either the Wulf's electroscope or the simple electroscope, and "ladle" charge to the interior of the pail.

The charges should be taken from an anode battery (522 34) or a power unit (522 35 a). If the process is continued for a long enough period very high voltages can be attained.

### 3. Experiments

in the field of electrostatics (see DK-Section 537.2).

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