# Bandgenerator

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
Fig. 1


Fig. 2

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#### 1. Description

The apparatus described below is an electrostatic high-tension generator developed by van de Graaf. It operates on the principle of self-excitation, and when running idle it produces a D. C. voltage of 120,000 to 140,000 volts. At approximately  $10^{-6}$  Coulomb, the loads are so small that the generator can be touched without involving any risk. When short-circuited, it produces a current of approximately  $10 \mu\text{A}$ .

The van de Graaf generator consists of a hollow metal sphere 18 cm in diameter, mounted on two plexiglass insulators which are about 30 cm high. The plexiglass mountings are fitted in a sturdy metal base. A hole is bored in the bottom of the metal sphere through which a 5 cm wide, endless rubber band is passed. The band runs over a plastics roller and three metal rollers. The bottom metal roller can be operated by means of a gear transmission turned with a crank. Several metal shearings are used to

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#### Notes

1. The five-figure numbers quoted in brackets refer to the catalogue numbers of the respective apparatus.
2. The data DC . . . refer to the experiment descriptions published as a collection of leaflets in the "LEYBOLD PHYSIKALISCHE HANDBLÄTTER".
3. The data and illustrations are not binding in every detail for the design of the apparatuses. It is our sole aim always to adapt our manufacturing programme to the most recent knowledge gained in all scientific and technical fields.

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give the band the required load. The large sparking ball of the generator consists of two sections. The top section can be removed to give an explanation of the apparatus and to exchange the band.

Two 4 mm sockets for taking off voltage are fitted, one on the ball and the other on the base. It is commendable to earth the lower socket, base and crank.

## 2. Operation and Circuit (Fig. 2)

The rubber band generally takes on a negative charge as the result of rubbing against the plastics exciter roller (1). It transports this charge to the interior of the ball (K) where it is completely taken up by the metal rollers (3 and 4), and fed to the surface of the ball. On turning the crank, the ball charge continues to increase until a limit value is reached determined by losses due to leakage and brush-discharge. Occasionally the band and ball may also take on a positive charge, particularly at high relative humidity.

The charging process is accelerated by means of induction effected by locating at suitable points 3 shearings (a, b and c) with their sharp edges facing the rubber band. The positive roller (1) induces a negative charge in the shearing (a), which increases the band charge as a result of brush-discharge. A positive charge is induced on shearings b and c by the upward-running negative band. This charge is brush-discharged on the downward-running band. Consequently the negative charging of the ball and the positive charge of the exciter roller are further increased.

To accelerate the charging process and to increase the current take-off, a large charging density of the rubber band is needed. This is obtained by making the two sections of the band run close to each other between the two guide rollers (1 and 4).

## 3. Starting up

Provided the air is dry, the modified van de Graaf generator mostly starts up immediately when the crank is rotated uniformly. Greater resistance to rotation is felt as the result of increased friction caused by the electrostatic attraction of the two sections of the band.

To make sure that the generator will start to operate, it is advisable to

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adjust the shearing (b) so that it is set at an angle and pointing upwards. This shearing can be made to rotate around its point of location. It is only when a very high voltage is required that the shearing should be directed downwards at an angle because then the conductance of a disruption spark along the band from the ball to earth is made difficult.

If the generator does not start up at once or only supplies low voltage, this condition can be remedied as follows: During rotation gently touch the exciter roller (1) with the dry hand. As a result the exciter roller is charged. Another method is to slightly press the rubber band with a dry cloth against the shearing (b).

If the generator still fails to supply voltage, then the rubber band or the plexiglass insulators are in all probability not completely dry. This is often the case in wet weather or after the generator has been carried from a cold storage room into the warm lecture room. The moisture can be removed by warming the apparatus with some form of heater (radiator or hot-air blower).

Instead of being operated by hand, the generator can also be driven by means of the rotator (Cat. No. 347 25 a). For this purpose, remove the nut on the crank, and screw on in its place the flexible coupling (Cat. No. 541 75). The other end of the coupling should be clamped in the gripper socket of the rotator, set at the right height.

#### 4. Maintenance

The high-tension ball, the rubber band and the plexiglass insulators must be kept perfectly clean otherwise the obtainable voltage will drop considerably. The band may be cleaned with acetone (free from water); powdering is harmful.

After a long period of operation, it may prove advisable to change the band round so that the under surface is facing upwards. To do this, loosen the lock nuts, and temporarily remove either roller (2) or roller (3). Spare bands are also supplied when needed.

The rollers (2) and (3) are somewhat thicker in the centre than at the edges to ensure that the band runs straight. However, if the band runs off the rollers after a long period of service or following an exchange of the band, then the roller axles can be re-centred because rollers (1), (2) and (3) run on short pins which are located

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eccentrically on the mounting pivots. After loosening the lock nuts, the position of the rollers can be adjusted by means of a screwdriver. The correct position will be found after a few trials, and then the lock nuts must be tightened up again.

## 5. Experiments

Electric field around a sphere  
DC 537.211

Electric potential lines around a sphere DC 537.213

Attraction of a charge by an earthed sphere (Image force) DC 537.217;b

Adjustment of image force in the electric field of a sphere DC 537.217;c

Measuring the amount of electricity in a charged sphere DC 537.214

In addition to these experiments, the van de Graaf generator may be used as an auxiliary apparatus in most experiments on electrostatics, e. g.

Charging of high-tension condensers (538 82) and Leyden jars (544 02)

Source of voltage for most of the experiments carried out with Wulf's electroscope (546 01 a)

Operation of discharge tubes and small X-ray tubes

Source of tension for the field emission electron microscope (554 60).

The way in which the van de Graaf generator operates can be demonstrated extremely well from the creepage sparks along the band, and with the aid of the luminous and brush-discharge phenomena visible in a completely blacked-out room.

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