P73.7

Seignettesalz-Kristall in Kapsel

Rochelle Salt Crystal in Case

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The Rochelle Salt Crystal is used to demonstrate piezoelectric and ferroelectric phenomena and to make use of them for measuring purposes.

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1. Description

In the middle of a transparent plastic box approx. 70 mm x 60 mm x 30 mm in size a quadratic plate of the crystal of approx. 4 mm thickness is arranged edgewise between the top and bottom plate of the box. The crystal consists of potassium-sodium tartrate (KOOOC-CHOH-CHOH-COONa · 4 H₂O). The two greatest areas of the crystal are about 7 cm² in size each and are covered with aluminium foil, and connected with the two 4 mm insulated sockets on the sides of the box by means of tin foil strips. In order to damp mechanical effects, the top and bottom outsides of the box are provided with a rubber layer.

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Due to its fragility the crystal should not be exposed to loads exceeding 5 kg.wt.

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2. Principe

By pressing or extending the crystal in the direction of one polar axis the atoms are displaced within the unit cells, and the mutual compensation

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Notes

- The five-digit numbers quoted in parentheses refer to the catalogue numbers of the respective apparatus.
- The specifications and illustrations are not binding in every detail for the design of the apparatus. It is our policy always to keep our manufacturing programme right up to date so that it makes full allowance for the latest knowledge acquired in all scientific and technical fields.

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of the atomic charges is interfered with, so that charges are set free on the surface. This phenomenon, which is called "piezoelectric effect", is exhibited by many crystals, and very pronouncedly so in case of Rochelle salt crystals. So this substance is particularly suitable for the investigations. On the two sides of the crystal, at right angles to the direction of pressure, charges of opposite sign are produced. The magnitude of these charges depends upon the size of the crystal surfaces and upon the pressure, but is independent of the thickness of the crystal.

In the direction of the crystal axis, the dielectric susceptibility is extremely high. On this account a very great polarization of the crystal can be brought about by applying an external electric field. This phenomenon is called "ferroelectric", for when the field is switched off, a remanent dielectric polarization of the crystal persists — in analogy to the remanent magnetism of a ferromagnetic substance.

3. Use

The crystal used here, when exposed to a force of 5 kg.wt., will carry a surface charge density of approx. $1,2 \times 10^{-10}$ As/cm². The low voltage which is thereby set up across the terminals of the box can be demonstrated, for example, with the Wulf electroscopes (546 00), or with the measuring amplifier (532 01). If the crystal is exposed to a varying pressure, using, for example, a sounding tuning fork (411 81 or 414 42), an oscillating voltage will be obtained, which can be demonstrated with the school cathode ray oscilloscope (575 15) or with the amplifier (580 00) with a speaker (580 05 or 580 06) in circuit.

In most experiments the box with the crystal is placed on the table and weighted with a weight or a tuning fork, or compressed by hand against its rubber plates. Horizontally directed impulses are best applied with the use of the prism table (460 25).

4. Experiments

a) Demonstration of the piezoelectric effect:

When the counter-electrode is brought into a position where the electroscopes is very sensitive, the leaf of the Wulf electroscopes will start to twitch as soon as the crystal connected to it is subjected to a sudden pressure pulse.

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With the measuring amplifier, range 3×10^{-7} As, one can show in addition that the charges produced are proportional to the changes of pressure. For this put on a 200 g weight, for example, at varying speeds and then put on weights from 50 to 1000 g one after the other and evenly, and observe the deflection. The latter experiment will become even more accurate on observing the deflection when removing the weights (this can be done more evenly). (See the book "Meßverstärker", Cat.-No. 532 04, expt. 27).

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The wave forms produced under the most various experimental conditions can be observed with the cathode ray oscilloscope (see the booklet: Ernst, "Versuche mit dem Kathodenstrahl-Oszillographen", Cat.-No. 575 16, expt. 131).

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b) When measuring the velocity of sound in solids with the demonstration scaler (559 62), the sound pulse produced by a bifilarly suspended steel ball impinging against the end of a metal rod is transformed into an electric pulse by the Rochelle salt crystal arranged at the other end of the rod (see Fig. 2 and the Directions for Use of the demonstration scaler 559 62).

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
c) Demonstration of the ferroelectric effect:

The school cathode ray oscilloscope can be used to show a hysteresis loop. For this make connections as shown in Fig. 3. The mode of action of this circuit becomes clear if the $0,2 \mu\text{F}$ capacitor is not considered initially. In that case the fraction of the mains A.C. derived from the 4200Ω potentiometer (537 01) is supplied to the pair of plates that effect the horizontal deflection of the beam. The full A.C. voltage of the mains is maintained across the crystal. The polarization

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current flowing in this circuit is supplied to the $2 \mu\text{F}$ capacitor where it will build up a voltage that is proportional to the polarization of the crystal. Through the amplifier of the oscilloscope, this voltage is supplied to the pair of plates that effect the vertical deflection of the beam. In this manner the polarization is represented as a function of the polarizing voltage.

If the hysteresis loop becomes too small, the $2 \mu\text{F}$ capacitor should be exchanged for a $1 \mu\text{F}$ one, and if it becomes too great, a $3 \mu\text{F}$ capacitor should be used instead.

The $0,2 \mu\text{F}$ capacitor has the sole function of compensating an undesirable phase shift that would otherwise occur in this circuit.

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