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(Federwaagen)

Prof. Maey's Dynamometers

(Spring Balances)

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(Balances à ressort)


Fig. 1 (1 : 10)

Kraftmesser 314.11 und 314.20 an Stativmaterial (2 Stativfüße 300.02, 2 Stativstangen 300.42, 2 LEYBOLD-Muffen 301.01)

Dynamometers 314.11 and 314.20 on stands (2 stand bases 300.02, 2 stand columns 300.42, 2 LEYBOLD's universal clamps 301.01)

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1. Description

Prof. Maey's dynamometers are spring balances, equipped at one end with a suspension ring and at the other with a hook. The springs are enclosed in a cardboard or light-metal sleeve and are surrounded by an easily visible, red and white scale. When the springs are weighted they expand, and the scale slides out of the sleeve thus enabling a reading to be made from a distance.

The dynamometers are supplied with a measuring range of 1 p; 2 p; 5 p; 10 p; 100 p; 1 kp; 2 kp; and 10 kp ($1 \text{ kp} = 1 \text{ kilopond} = 1000 \text{ p} = 9,81 \text{ kg} \cdot \text{m/sec}^2 = 9,81 \text{ Newton}$).

2. Usage

Although the springs are extremely durable, they should not be loaded with weights heavier than those for which they are designed, in other

Notes

1. The five-figure numbers quoted in brackets refer to the catalogue numbers of the respective apparatus.
2. The data DK . . . refer to the experiment descriptions published as a collection of leaflets in the "LEYBOLD PHYSIKALISCHE HANDBLÄTTER".
3. The data and illustrations are not binding in every detail for the design of the apparatuses. It is our sole aim always to adapt our manufacturing programme to the most recent knowledge gained in all scientific and technical fields.

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words, so that they do not expand beyond the end of the scale. Moreover the springs should be held by the ring only, and not by the sleeve which is not suited for this purpose but merely designed as a guide.

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3. Experiments

(see DK-Sections 531.2 to 531.6)

The dynamometers enable the explanation of the term force, the demonstration of composition and resolution of forces as well as the quantitative handling of such simple machines as inclined plane (341 21), lever (342 58), and wheel and axle (342 75). Furthermore they may also be used as provisional balances.

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